Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room W-066-
Silver Spring, MD 20993-0002

SMT Schilling Metalltechnik GmbH
% Medagent GmbH & Co. KG
Mr. Erik Schilling
Regulatory Affairs Manager
Griesweg 47
Müehlheim, Baden-Württemberg
Germany 78570

SEP 10 2010

Re: K100736

Trade/Device Name: SMT Schilling Kirschner/Guide Wires
Regulation Number: 21 CFR 888.3040
Regulation Name: Smooth or threaded metallic bone fixation fastener
Regulatory Class: II
Product Code: HTY, JDW
Dated: September 2, 2010
Received: September 7, 2010

Dear Mr. Schilling:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21

CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Mark N. Melkerson
Director
Division of Surgical, Orthopedic
and Restorative Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

K100736

510(k) Number: K100736

Device Name:
SMT Schilling Kirschner/Guide Wires

SEP 10 2010

Indications For Use:

Orthopaedic fixation pins and wires are intended to perform as fixation and stabilization unit of bone fractures or as guidance at insertion of implants into the skeletal system.

Prescription Use YES
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Jonita J. for mxm
(Division Sign-Off)
Division of Surgical, Orthopedic,
and Restorative Devices

Page 1 of 1

510(k) Number K100736

K100736

510(K) NOTIFICATION

PREMARKET NOTIFICATION [510(K)] SUBMISSION

DATE OF APPLICATION: 2010-07-16

APPLICANT: SMT Schilling Metalltechnik GmbH
Griesweg 33
78570 Mühlheim an der Donau
Germany
Tel.: (07463) 99309-0
Fax: (07463) 99309-59

SIGNATURES:

Erik Schilling, Managing Director

Manfred Huber, Quality Manager

Table of Contents

1. Establishment Registration Number 3
 2. FDA-Designated Agent..... 3
 3. Device Name 3
 4. Description of the Device..... 3
 5. Intended Use 3
 6. Classification Product Code / Subsequent Code 4
 7. Performance Standards..... 4
 8. Labeling / packaging..... 5
 8.1. Sterile products..... 5
 8.2. Shelf-Life of sterile Packed Devices 5
 8.3. Non sterile products..... 6
 9. Biocompatibility 6
 10. Sterility 7
 10.1. Non sterile Products..... 7
 10.2. Sterile Products..... 7
 11. Certified Quality Management System 7
 12. Substantial Equivalence Comparison..... 8

Appendices

Annex	Description of Annex	Total Pages
I	Truthful and Accurate Statement	1
II	510(k) Summary of Safety and Effectiveness in Accordance with SMDA 1990	2
III	IFU-Statement	1
IV	Instructions for Use	6
V	Catalogues	6
VI	Test Reports	16
VII	Predicate Devices; Instruction for Use, 510(k) Summaries, Catalogue, Literature	20
VIII	Engineering Drawings	12
IX	Standards Forms (Form 3654)	

1. Establishment Registration Number

SMT Schilling Metalltechnik GmbH: 9680619

2. FDA-Designated Agent

Sanjiv Kumar
Myco Medical Supplies Inc.
158 Towerview Court
Cary, NC 27513
U.S.A.
Phone: 919 4602535 Ext
Fax: 919 4602536
Email: Skumar@Mycomedical.Com

3. Device Name

Trade Name: Kirschner Wire / Guide Wire
Common Name: Kirschner Wire (K Wire)
Classification Name: Pin, Fixation

4. Description of the Device

Orthopaedic fixation pins and wires are metal pins for use in fixation of bone fractures, for bone reconstructions, as guide pins for insertion of other implants or implantation through the skin so that traction may be applied to the skeleton system.

To ensure the multi-use of these devices, many different models are available. The differences can be as follows:

- Diameter: from 0.6 up to 6.35mm (0.020 up to 0.250 inch)
- Length: from 60 up to 500 mm (2.36 up to 19.69 inch)
- Tips: diamond or trocar Point, round, flat, with or without 3- or 4- shank ends, with or without spherical shape.
- Surface: complete or partial smooth and / or threaded, with or without threading cutter.

5. Intended Use

Orthopaedic fixation pins and wires are intended to perform as fixation and stabilization unit of bone fractures or as guidance at insertion of implants into the skeletal system.

6. Classification Product Code / Subsequent Code

Orthopaedic fixation pins and wires can be classified according following Device Names and Product Codes:

Device:	Pin, Fixation Smooth	Pin, Fixation, Threaded
Medical Specialty:	Part 888, Orthopedic	Part 888, Orthopedic
Product Code:	87 HTY	87 JDW
Device Class:	2	2
Regulation Number:	888.3040	888.3040

Description according 21CFR888.3040:

"Sec. 888.3040 Smooth or threaded metallic bone fixation fastener.

(a) Identification. A smooth or threaded metallic bone fixation fastener is a device intended to be implanted that consists of a stiff wire segment or rod made of alloys, such as cobalt-chromium-molybdenum and stainless steel, and that may be smooth on the outside, fully or partially threaded, straight or U-shaped; and may be either blunt pointed, sharp pointed, or have a formed, slotted head on the end. It may be used for fixation of bone fractures, for bone reconstructions, as a guide pin for insertion of other implants, or it may be implanted through the skin so that a pulling force (traction) may be applied to the skeletal system.

(b) Classification. Class II."

7. Performance Standards

Orthopaedic fixation pins and wires are made out of following raw materials, which are meeting both American and International Standards:


Implant Grade Stainless Steel

Meets all requirements of the special steel grade 316 L according to ASTM F 138 / ASTM F 139 as well as ISO 5832-1.


8. Labeling / packaging

8.1. Sterile products

Sterile products are primarily packed in plastic blisters and Tyvek lids, secondarily packed in PU plastic bags likewise follows:

Blister packed in cardboard box with attached label:

Each delivered pack is labeled as follows (on cardboard box):

REF 31-61-23-1635		K-Wire		Material: 1.4441 - Stainless Steel		CE 043		Rx ONLY		STERILE R		QTY 1		V1488 10060	
SMT Schilling Metalltechnik GmbH Grossegang 33 D-78570 Mühlheim a.d. Donau		2010-03 Manufacturing date		2015-02 Best before		SMT SCHILLING GmbH Metalltechnik		for single use only		sterile, sealed		sterilizing documents		REF 31-61-23-1635	
K-Wire		Material: 1.4441 - Stainless Steel		STERILE R		QTY 1		REF 31-61-23-1635		K-Wire		Material: 1.4441 - Stainless Steel			
REF 31-61-23-1635		K-Wire		Material: 1.4441 - Stainless Steel		SMT SCHILLING GmbH Metalltechnik		V1488 10060		STERILE R		QTY 1			
SMT Schilling Metalltechnik GmbH Grossegang 33 D-78570 Mühlheim a.d. Donau		2010-03 Manufacturing date		2015-02 Best before		SMT SCHILLING GmbH Metalltechnik		for single use only		sterile, sealed		sterilizing documents		LOT 1009001	
K-Wire		Material: 1.4441 - Stainless Steel		STERILE R		QTY 1		REF 31-61-23-1635		K-Wire		Material: 1.4441 - Stainless Steel		Rx ONLY	

8.2. Shelf-Life of sterile Packed Devices

Due to an accelerated ageing simulation according to ASTM F 1980-02, a shelf life of 5 Years can be determined. See test report related accelerated aging testing in appendix A06.

8.3. Non sterile products

Non sterile products are packed in tubes made of transparent, rigid PVC. One tube can contain up to 10 pieces according to customers needs.

9. Biocompatibility

The biocompatibility is guaranteed by the composition of the used materials as mentioned in chapter 6. These materials are also used for most other surgical implants as bone screws, bone plates or neuro surgical implants as Aneurysm Clips as mentioned in following 510(k)'s:

- K983758, Aesculap AG, Aesculap Yasargil Titanium Aneurysm Clips;
- K000080, Howmedica Osteonics Corp., Ansis III Cannulated Screw System;
- K003500, Rebstock Instruments GmbH, Rebstock Yasargil Aneurysm Clips.

SMT Schilling Metalltechnik GmbH's Bone Wires are only made out of the in chapter 6 mentioned Raw Materials.

Compared with competitors SMT Schilling Metalltechnik GmbH's Bone Wires are made out of similar materials.

10. Sterility

The devices are delivered in either sterile or non-sterile conditions.

10.1. Non sterile Products

The user may sterilize these devices by using a validated and applicable sterilization process.

SMT Schilling Metalltechnik GmbH recommends using a steam-sterilizer that uses a validated sterilization cycle of 134°C / 273°F, pre-vacuum, for 4 minutes and 15 minutes minimum drying time.

The steam sterilization process will be validated according to ISO 17665-1 with a SAL of 10^{-6} .

10.2. Sterile Products

The gamma radiation sterilization process for sterile products has been validated under consideration of an SAL 10^{-6} according to the requirements based on ISO 11137 (see appendix A06) using an irradiation dose of 5.0 kGy to ensure effective conditions ($SAL \leq 10^{-6}$) for the sterilization of the products by application of a minimum dose of 25 kGy.

The tested products on the test report **082248-10; TPLO Plate left and right** have been identified as worst case products based upon dimensions and geometry. Therefore an effective sterilization method can be determined for Kirschner / Guide Wires. The tested material (stainless steel 1.4441) is referenced on the drawing of respective article no. attached to the test reports.

The packaging of sterile products have been validated under consideration of ISO 11607-1 and 2 (see appendix A06).

11. Certified Quality Management System

SMT Schilling Metalltechnik GmbH is certified according following Standards or Normative Directives:

- DIN EN ISO 13485:2003
- European Medical Device Directive (MDD) 93/42/EEC Annex II.

The certified system is continuously checked by Notified Body MDC GmbH, CE-Code # 0487

12. Substantial Equivalence Comparison

Compared with competitors Bone Wires, SMT Schilling Metalltechnik GmbH's Bone Wires are similar / Different in following areas:

Area / Competitor	SMT Metalltechnik GmbH Orthopaedic Fixation Pins and Wires	Internal / external fixation devices (#K070561)	Störk Kirschner Wires (K-Wires) and Steinmann Pins (#K030665)
indication for use	Orthopaedic fixation pins and wires are intended to perform as fixation and stabilization unit of bone fractures or as guidance at insertion of implants into the skeletal system.	The KMedic® Internal/External Fixation Devices are non-sterile, single-use, unilateral internal/external internal/external fixation devices intended to be used for fixation in the treatment of bone conditions including limb lengthening, osteotomies, arthrodesis, fracture fixation and other bone conditions amenable to treatment by use of the internal/external fixation modality. similar	The intended use of Störk Instrumente GmbH's Bone Wires is the fixation of bone fractures and / or to guide other implants during insertion to the skeletal system. Due to their design and used materials, which is either Stainless Steel acc. ASTM F 138 / 139 (316L) or Titanium Alloy acc. to ASTM F 136, Grade 5, these devices can be used as implants. Störk Instrumente GmbH is offering their bone wires only in Non-Sterile condition. similar
design	Pins with smooth or threaded surface and various types of tips and ends	Pins and Wires of various lengths, sizes and end configurations similar	Pins with smooth or threaded surface and various types of tips and ends similar
materials / biocompatibility	ASTM F 138 / 139 (316L)	ASTM F 138 / 139 (316L) Equivalent	ASTM F 138 / 139 (316L) or Titanium Alloy acc. to ASTM F 136, Grade 5 Equivalent
sterility	Delivered in sterile or non-sterile conditions	Delivered in non-sterile conditions. similar / DIFFERENT ¹	Delivered in non-sterile conditions. similar / DIFFERENT ¹
Dimensions	Diameter: ø0,6 up to ø6,35 mm Length: 60 up to 500 mm	Diameter: ø0,7 up to ø4,5 mm Length: 101,6 up to 304,8 mm Similar	Diameter: ø0,5 up to ø6,5 mm Length: 60 up to 500 mm Similar

¹SMT Schilling's devices are delivered in either sterile or non-sterile conditions. For verification of sterilization and sterile packaging, see annex VI.

Truthful and Accurate Statement

**PREMARKET NOTIFICATION
TRUTHFUL AND ACCURATE STATEMENT**
[As required by 21 CFR 807.87(j)]

I certify that, in my capacity as a **Managing Director** of **SMT Schilling Metalltechnik GmbH**, I believe to the best of my knowledge, that all data and information submitted in the premarket notification are truthful and accurate and that no material fact has been omitted.

(Signature)

Erik Schilling, CEO

2010/07/16
(Date)

K100736
(Premarket Notification [510(k)] Number)

K100736

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS IN ACCORDANCE WITH SMDA OF 1990

Date of Application: 2010-07-16

APPLICANT: SMT Schilling Metalltechnik GmbH
Griesweg 33
78570 Mühlheim an der Donau
Germany
Tel.: (07463) 99309-0
Fax: (07463) 99309-59

CONTACT PERSON: Mr. Erik Schilling
Managing Director
Tel.: (07463) 99309-0
Fax: (07463) 99309-59

1. Device Name

Trade Name: Orthopaedic Fixation Pins and Wires / Kirschner / Guide Wires
Common Name: Kirschner Wire (K-Wire)

2. Classification Product Code / Subsequent Code

Our implant system can be classified according following device names and product codes:

Device:	Pin, Fixation Smooth	Pin, Fixation, Threaded
Medical Specialty:	Part 888, Orthopedic	Part 888, Orthopedic
Product Code:	87 HTY	87 JDW
Device Class:	2	2
Regulation Number:	888.3040	888.3040

3. Substantial Equivalence

Orthopaedic fixation pins and wires are substantial equivalent based upon design, dimensional and materials characterization to the Störk Kirschner Wires (K-Wires) and Steinmann Pins (#K030665) of Stork Instrumente GmbH, 78576 Emmingen-Liptingen, Germany and Teleflex KMedic Internal/External Fixation Devices (#K070561) of Teleflex Medical, Bannockburn, IL.

4. Description of the Device

Orthopaedic fixation pins and wires are metal pins for use in fixation of bone fractures, for bone reconstructions, as guide pins for insertion of other implants or implantation through the skin so that traction may be applied to the skeleton system.

To ensure the multi-use of these devices, many different models are available. The differences can be as follows:

- Diameter: from 0.6 up to 6.35mm (0.020 up to 0.250 inch)
- Length: from 60 up to 500 mm (2.36 up to 19.69 inch)
- Tips: diamond or trocar Point, round, flat, with or without 3- or 4- shank ends, with or without spherical shape.
- Surface: complete or partial smooth and / or threaded, with or without threading cutter.

5. Intended Use

SMT Schilling's orthopaedic fixation pins and wires are intended to perform as fixation and stabilization unit of bone fractures or as guidance at insertion of implants into the skeletal system.

6. Conclusion

Based on the available 510(k) summaries and 510(k) statements (21 CFR 807) and the information provided herein, we conclude that SMT Schilling's Kirschner Wires and Pins are substantially equivalent to the existing legally marketed devices under Federal Food, Drug and Cosmetic Act.

Indications for Use

510(k) Number: **K100736**

Device Name:

SMT Schilling Kirschner / Guide Wires

Indications for Use:

Orthopaedic fixation pins and wires are intended to perform as fixation and stabilization unit of bone fractures or as guidance at insertion of implants into the skeletal system.

Prescription Use **YES**
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use **NO**
(21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

SMT Metalltechnik GmbH
Schilling Orthopedic Surgical Instruments
Griesweg 33
D – 78570 Mühlheim/Donau
Fon: 07463 – 99309-0
E-Mail: info@smt-metalltechnik.de
www.smt-metalltechnik.de

INSTRUCTIONS FOR USE (IFU – 011.1)

English

Instructions for Use: Orthopaedic Fixation Pins and Wires (Guide Wires / Kirschner Wires)

Contents

The package contains one implant

Material

The implant is made by following materials:

- Stainless steel (ASTM F 138 / ASTM F 139 as well as ISO 5832-1).

It is supplied either sterile or non sterile and available in numerous sizes

Intended use

Indications

Orthopaedic fixation pins and wires are intended to perform as fixation and stabilization unit of bone fractures or as guidance at insertion of implants into the skeletal system.

SMT Metalltechnik GmbH
Schilling Orthopedic Surgical Instruments
Griesweg 33
D – 78570 Mühlheim/Donau
Fon: 07463 – 99309-0
E-Mail: info@smt-metalltechnik.de
www.smt-metalltechnik.de

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Contraindications

- Do not use the implant in cases of:
- States of health, which exclude adequate assistance or healing process, i.e.:
 - Impairment of Blood supply
 - Inadequate quantity / quality of bone structure
 - Extreme adiposis
 - previous infection
 - Strong twist or disposition of shank
 - Mental conditions which exclude participation in rehabilitation programs (Parkinson's, alcohol or drug abuse etc.)
 - Major physical activities, adherent with intense percussion, on which the implants are exposed to excessive pressure.
 - Allergic reactions against components.

The pins and wires have not been evaluated for safety and compatibility in the MR environment. The pins and wires have not been tested for heating or migration in the MR environment. Therefore do not use the pins and wires in a MR environment.

Adverse reactions

- Adverse reactions may include:
- Clinical failure (i.e. pain or injury) due to bending, loosening, wear and tear, fracture of implant, loss of fixation, dislocation and/or migration
 - Pain and/or abnormal sensations due to the presence of the implant
 - Primary and/or secondary infections
 - Allergic reactions to implant material
 - Injury to vessels, nerves and organ
 - Hematoma and /or impaired wound healing; hemorrhage

Safety precautions

- Prior to use, thoroughly read these instructions for use and become familiar with the surgical technique.
- Keep the instructions for use accessible to all staff.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques. Proper surgical performance of the implantation is the responsibility of the operating surgeon.
- The manufacturer is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrect operating techniques, the limitations of treatment methods or inadequate asepsis.

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- The operating procedure must be explained to the patient, and the patient's understanding of the following information must be documented:
- The patient is aware of the risks associated with general surgery, orthopedic surgery and with general anesthesia.
- The patient has been informed about the advantages of the implant procedure and about possible alternative treatments.
- The implant can fail due to excessive load, wear and tear, or infection.
- The service life of the implant is determined by body weight and physical activity. The implant must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- Corrective surgery may be necessary if the implant fails.
- The patient must have their physician carry out follow-up examinations of the implant at regular intervals.

Intraoperative

- Prior to use, verify the integrity of the sterile packaging and check the product expiration date.

WARNING:

Never use implants if the packaging is damaged. Never use implants that are past their expiration date without resterilizing.

Postoperative

- Reiterate preoperative instructions to the patient.
- Ensure that the patient is aware of physical activity restrictions and possible adverse reactions.

Warranty

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with the instructions for use.
 Technical alterations reserved.

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


Wires-Article No z.B 31-61-21-2201
Length of the wires as 21cm=21
Diameter as 2,2mm=22
Style as-01=Trokar-flat

Trokar Standard Winkel w 15° bis Ø4,8
20° bei Ø6,35
Trokar standard angle 15° to up to Ø4,8
20° at Ø6,35

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Drawing# = item#
 Direct Reference to
 Catalogue

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SMT Schilling Metaltechnik GmbH