SMT Metalltechnik GmbH
Schilling Orthopedic Surgical Instruments
Griesweg 33
D – 78570 Mühlheim/Donau
Fon: 07463 – 99309-0
E-Mail: info@smt-metalltechnik.de
www.smt-metalltechnik.de

INSTRUCTIONS FOR USE (IFU – 011.1)

English

Instructions for Use: Orthopaedic Fixation Pins and Wires (Guide Wires / Kirschner Wires)

Contents

The package contains one implant

Material

The implant is made by following materials:

- Stainless steel (ASTM F 138 / ASTM F 139 as well as ISO 5832-1).

It is supplied either sterile or non sterile and available in numerous sizes

Intended use

Indications

Orthopaedic fixation pins and wires are intended to perform as fixation and stabilization unit of bone fractures or as guidance at insertion of implants into the skeletal system.

SMT Metalltechnik GmbH
Schilling Orthopedic Surgical Instruments
Griesweg 33
D – 78570 Mühlheim/Donau
Fon: 07463 – 99309-0
E-Mail: info@smt-metalltechnik.de
www.smt-metalltechnik.de

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Contraindications

Do not use the implant in cases of:

- States of health, which exclude adequate assistance or healing process, i.e.:
 - Impairment of Blood supply
 - Inadequate quantity / quality of bone structure
 - Extreme adiposis
 - previous infection
 - Strong twist or disposition of shank
- Mental conditions which exclude participation in rehabilitation programs (Parkinson's, alcohol or drug abuse etc.)
- Major physical activities, adherent with intense percussion, on which the implants are exposed to excessive pressure.
- Allergic reactions against components.

The pins and wires have not been evaluated for safety and compatibility in the MR environment. The pins and wires have not been tested for heating or migration in the MR environment. Therefore do not use the pins and wires in a MR environment.

Adverse reactions

Adverse reactions may include:

- Clinical failure (i.e. pain or injury) due to bending, loosening, wear and tear, fracture of implant, loss of fixation, dislocation and/or migration
- Pain and/or abnormal sensations due to the presence of the implant
- Primary and/or secondary infections
- Allergic reactions to implant material
- Injury to vessels, nerves and organ
- Hematoma and /or impaired wound healing; hemorrhage

Safety precautions

- Prior to use, thoroughly read these instructions for use and become familiar with the surgical technique.
- Keep the instructions for use accessible to all staff.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques. Proper surgical performance of the implantation is the responsibility of the operating surgeon.
- The manufacturer is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrect operating techniques, the limitations of treatment methods or inadequate asepsis.

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- Under no circumstances may modular implant components from different suppliers be combined.
- Each patient's record shall document the implant used (name, article number, lot number).
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed about post- surgical behavioral requirements.
- Damage to the weight-bearing structures can give rise to loosening, dislocation and migration, as well as other complications. To ensure the earliest possible detection of such catalysts of implant dysfunction the implant must be checked periodically post operative using appropriate techniques
- Never reuse an implant. Although the implant may appear undamaged, previous stresses may have created non-visible damage that could result in implant failure
- Never use implants if the packaging is damaged. An implant with damaged packaging might be damaged and thus may not be resterilized/used.
- Never use implants that are past their expiration date without resterilizing.
- Do not resterilize implants if they have been contaminated by blood or secretions.

Storage, inspection and sterilization

Storage

The implant is individually packed in protective packaging that is labeled according to its contents.

The implant is sterilized with gamma sterilization (25kGy minimum)

- Always store the implant in the original protective packaging.
- Do not remove the implant from the packaging until immediately before use.
- Store the implant in normal hospital environmental conditions.

Cleaning / Disinfection

Products delivered in non-sterile condition or sterile products, which sterile packaging has been opened or damaged, must be cleaned, disinfected, and sterilized prior to use. For cleaning and sterilization remove the implant from its packaging since the packaging may not be used to clean, disinfect and sterilize the product. A suitable cleaning, disinfection and sterilization process must be applied by the user. Only pH-neutral cleaning agents should be used. The preparation references of the respective cleaning and disinfection agent manufacturer must be considered.

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Sterilization / Resterilization

Non sterile implants as well as Implants from open packages which have not been implanted may be sterilized by the user. The implants must be sterilized using a process that has been validated by the health care provider – reference ISO 11134, “Sterilization of health care products – Requirements for validation and routine control – Industrial moist heat sterilization” and AAMI TIR 12, “Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: a guide for device manufacturers”. SMT has shown that the implant can be sterilized using the following steam sterilization cycle:

- Pre-vacuum steam sterilization, wrapped items
- 134 °C for 4 minutes, minimum Dry Time 15 minutes

The health care provider must validate that the steam sterilization cycle is effective in their facility since all steam sterilization chambers are unique.

Warning:

the device may only be resterilized if it has not been in touch with bodily fluids

Procedure

Preoperative

The operating surgeon draws up an operation plan specifying and documenting the following:

- Implant component(s) and their dimensions.
- Proper position of the implant component(s)

The following conditions must be fulfilled prior to application:

- All required implant component(s) are readily available.
- Highly aseptic operating conditions are present.
- All requisite implantation instruments must be available and in working order, including the trial

Warning:

Never use or process damaged or defective instruments. Contact your local representative or dealer for repair or replacement. The use of an instrument for tasks other than those for which they are intended may result in damaged / broken instruments or patient injury

- The operating surgeon and operating room team must be thoroughly familiar with the operating technique, as well as the range of implants and instruments to be applied. Complete information on these subjects must be ready available at the workplace.

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- The operating procedure must be explained to the patient, and the patient's understanding of the following information must be documented:
- The patient is aware of the risks associated with general surgery, orthopedic surgery and with general anesthesia.
- The patient has been informed about the advantages of the implant procedure and about possible alternative treatments.
- The implant can fail due to excessive load, wear and tear, or infection.
- The service life of the implant is determined by body weight and physical activity. The implant must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- Corrective surgery may be necessary if the implant fails.
- The patient must have their physician carry out follow-up examinations of the implant at regular intervals.

Intraoperative

- Prior to use, verify the integrity of the sterile packaging and check the product expiration date.

WARNING:

Never use implants if the packaging is damaged. Never use implants that are past their expiration date without resterilizing.

Postoperative

- Reiterate preoperative instructions to the patient.
- Ensure that the patient is aware of physical activity restrictions and possible adverse reactions.

Warranty

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with the instructions for use.

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Symbols / Labeling

The symbols used for labeling acc. to DIN EN 980, are described as follows:

	Single Use Only
	Attend instructions for use
	CE-Mark with number of notified body mdc Stuttgart
	Reference number
	Lot Number
	Non sterile
	Sterilized by radiation

For further information

Please contact SMT Schilling Metalltechnik GmbH or your authorized representative if further information on this product is needed

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Weitere Informationen erhalten Sie jederzeit bei SMT Schilling Metalltechnik GmbH oder Ihrem autorisierten Vertreter.
SMT Schilling GmbH; Griesweg 33; 78570 Mühlheim / Germany; www.smt-metalltechnik.de